

NirogScan Health Report

Patient Information

Name: asd
Age: 1
Gender: Male
Weight: N/A
Height: N/A
Notes: None


Recording & Vital Signs

ECG Duration: 33.2s
ECG Sample Rate: 125 Hz
Analysis Window: 22.82s - 27.82s
Window SNR: 21.3 dB
R-peaks Found: 5
Heart Rate (ECG): 58.1 +/- 1.7 bpm
SpO2: 81%
Pulse Rate (PPG): 51.7 bpm
Respiration Rate: 24.9 breaths/min
Perfusion Index: 204.52%


Generated: 2026-02-27 16:59:58

ECG Medical Grid Analysis

Cleaned ECG with R-peaks




Single Beat Detail




ECG PQRS Wave Detection


PQRST Points Labeled




Colored Segments



RR Tachogram



Poincare Plot



Heart Rate Variability (HRV) Analysis


Time Domain

Average RR Interval:	1034.00
Heart Rate Variability:	36.59
Short-term HRV:	47.33
Large Beat Changes %:	25.00
Small Beat Changes %:	25.00
Median RR Interval:	1020.00

Frequency Domain

Frequency analysis requires
60+ seconds of recording


RR Distribution



Nonlinear


Nonlinear analysis requires
10+ R-peaks

Full ECG Recording (33.2s)




PPG Signal Analysis


IR Channel




Red Channel



IR Baseline Tracking



Red Baseline Tracking



PPG Advanced Analysis

Oxygen Saturation

81.0%

SpO2


Age Index (Stiffness)

Age Index
not available


Pulse Arrival Time

PAT
not available


IR Instantaneous Frequency




Red Instantaneous Frequency



IR Beat Segmentation (5 beats)



Red Beat Segmentation (5 beats)



Respiration Rate

24.9

breaths/min

Baseline Tracking

IR Drift: -30.202/s
IR Std: 1904.75
Red Drift: -18.681/s
Red Std: 1874.22

Analysis Summary

ECG Analysis:

Heart Rate:	58.1 $\pm$ 1.7 bpm
PR Interval:	147 ms
QRS Duration:	90 ms
QTc:	390 ms
HRV_SDNN:	36.6
HRV_RMSSD:	47.3

PPG Analysis:

SpO2:	81%
Respiration:	24.9 br/min

Note: For informational purposes only. Consult healthcare professional for medical advice.

Understanding Your Results

(Simple explanations of each measurement)

Oxygen Saturation:

Blood oxygen level. Normal: 95-100%.

Perfusion Index:

Blood flow strength. Higher = better circulation.

Pulse Rate:

Heart rate from blood flow. Should match ECG heart rate.

Average RR Interval:

Average time between heartbeats. Normal: 600-1000ms.

Heart Rate Variability:

How much heartbeat timing varies. Higher = healthier heart. Normal: 50-100ms.

Short-term HRV:

Beat-to-beat changes. Reflects relaxation nervous system. Normal: 20-50ms.

Large Beat Changes %:

Percentage of beats differing >50ms. Higher = good health. Normal: 3-25%.

PR Interval:

Atrial to ventricular time. Normal: 120-200ms.

QRS Duration:

Ventricular contraction time. Normal: 80-120ms.

Corrected QT:

Heart rate adjusted QT. Normal: <440ms.

Note: These are for informational purposes only. Consult a healthcare professional for medical advice.