

NirogScan Comprehensive Health Report

Patient Information

Name: ak
Age: 24
Gender: Male
Weight: N/A
Height: N/A
Notes: None


Recording & Analysis Summary

ECG Duration: 11.8s (1475 samples)
ECG Sampling Rate: 125 Hz
Best Analysis Window: 1.98s - 6.98s
Window SNR: 20.6 dB
R-peaks Detected: 5
Analysis Engine: NeuroKit2 v0.2.12
ECG Heart Rate: 60.2 ± 2.1 bpm
PPG Duration: 23.6s


Report Generated: 2026-01-31 19:27:04

ECG Signal Analysis

Raw ECG (Analysis Window: 1.98s - 6.98s)




Cleaned ECG (nk.ecg_clean) - Quality: 0.484




ECG Morphological Analysis (PQRST Detection)


Single Beat with PQRST Waves




Average Beat Template (n=3)



RR Interval Tachogram



Poincaré Plot



HRV Analysis - Time Domain


Core Time Domain Metrics

Average RR Interval:	998.00
Heart Rate Variability:	38.30
Short-term HRV:	34.87
Percentage of Large Changes:	25.00
Percentage of Small Changes:	25.00
Successive Differences SD:	23.09


Additional Time Metrics

Median RR Interval:	1008.00
Median Absolute Deviation:	23.72
Coefficient of Variation:	0.04
Interquartile Range:	34.00
Triangular Index:	0.00
HRV Triangular Index:	4.00

RR Interval Distribution



Instantaneous Heart Rate



HRV Analysis - Frequency & Nonlinear Domains

Frequency Domain Metrics

HRV_TP: 0.00

Insufficient frequency data

Nonlinear Metrics (Poincaré)

Poincaré SD1: 16.3299

Poincaré SD2: 37.6652

SD1/SD2 Ratio: 0.4336

Poincaré Area: 1932.2994

Cardiac Sympathetic Index: 2.3065

Cardiac Vagal Index: 3.9930

Entropy & Fractal Metrics


Approximate Entropy: 0.4055

Sample Entropy: -inf

Fuzzy Entropy: 3.4816

PPG Signal Analysis (MAX30102 IR Channel)

PPG Signal with Peak Detection




Insufficient HR data

Insufficient data

Insufficient data

Full ECG Recording (11.8s total)



Analysis Summary

ECG Analysis:

Signal Quality (SNR):	20.6 dB
R-peaks Detected:	5
Heart Rate:	60.2 ± 2.1 bpm
SDNN:	38.3 ms
RMSSD:	34.9 ms

PPG Analysis:

Peaks Detected:	1
-----------------	---

Device Information:

ECG Sensor:	AD8232 (Gain: 1100)
ECG Sampling:	125 Hz
PPG Sensor:	MAX30102 (Red 660nm, IR 880nm)
PPG Sampling:	25 Hz
ADC:	12-bit, 3300mV ref

Understanding Your Results

(What each measurement means in simple terms)

- **Average RR Interval:**
The average time between heartbeats in milliseconds. Normal range: 600-1000ms.
- **Heart Rate Variability:**
How much your heartbeat timing varies. Higher is generally better, indicating a healthy, adaptable heart. Normal: 50-100ms.
- **Short-term HRV:**
Measures beat-to-beat changes. Reflects your parasympathetic (rest & digest) nervous system activity. Normal: 20-50ms.
- **Percentage of Large Changes:**
Percentage of heartbeats that differ by more than 50ms from the previous one. Higher values indicate good heart health. Normal: 3-25%.
- **Poincaré SD1:**
Short-term variability from Poincaré plot. Related to RMSSD.
- **Poincaré SD2:**
Long-term variability from Poincaré plot. Related to SDNN.
- **Approximate Entropy:**
Signal complexity/predictability. Lower = more regular heartbeat.
- **Sample Entropy:**
Similar to ApEn but more consistent. Measures heart rhythm complexity.
- **PR Interval:**
Time from atrial to ventricular activation. Normal: 120-200ms.
- **QRS Duration:**
Time for ventricles to contract. Normal: 80-120ms.
- **QT Interval:**
Total ventricular activity time. Varies with heart rate.
- **Corrected QT:**
QT interval adjusted for heart rate. Normal: <440ms.

*Note: These measurements provide insights into heart health but are not medical diagnoses.
Please consult a healthcare professional for medical advice.*